EmCyte Corporation
4331 Veronica S. Shoemaker Blvd.,
Suite 4 Fort Myers, Florida 33916
USA

PL-42-32
REV 8
Date: December 2018

AbsolutePRP™ 20mL
GenesisCS Component Concentrating System

PN: GS301-EU
Instruction for Use

ATTENTION OPERATING SURGEON

NOTE: DEVICE IS FOR SINGLE USE ONLY. Discard the entire disposable system after one use, using an acceptable disposal method for products potentially contaminated with blood.

DESCRIPTION

The Platelet Concentrating System 20mL is a concentrating system manufactured by EmCyte Corporation. The Platelet Concentrating System 20mL is part of the GenesisCS Component Concentrating System's family of kits. The kit concentrates whole blood components by density through the use of the 20mL Concentrating Device and the Executive Series Centrifuge II provided by EmCyte Corporation. The system contains syringes, a blood draw accessory and the Concentrating Device accessories. The Platelet Concentrating System 20mL permits the rapid preparation of platelet rich plasma from 20mL of the patient's whole blood drawn at the time of treatment.

MATERIALS

The materials used are syringes, needles, tubing, connectors, and concentrating devices. The materials consist of medical grade polymers, elastomers and stainless steel that are suitable for use in medical devices. All components in this system are packaged, labeled and sterilized as indicated by the manufacturer's labeling. All components in this system are latex-free.

INDICATIONS FOR USE STATEMENTS

FDA 510(k) BK050055

The GenesisCS Component Concentrating System is designed to be used for the safe and rapid preparation of autologous platelet rich plasma (PRP) from a small sample of blood at the patient's point of care. The PRP can be mixed with autograft and allograft bone prior to application to an orthopedic surgical site as deemed necessary by the clinical use requirements.

FDA 510(k) K07066

The GenesisCS Component Concentrating System is designed to be used in a clinical laboratory or intraoperatively at the point of care for the safe and rapid preparation of platelet poor plasma and platelet concentrate from a small sample of blood and for a preparation of a cell concentrate from bone marrow. The safety and effectiveness of this device for in vivo indications for use has not been established.

FDA 510(k) K024080

The SECQUIRE cell separator is designed for the safe and effective preparation of low volume platelet rich plasma, and platelet poor plasma at the point of care.

CE0653

The Concentrating System is designed to be used for the safe and rapid preparation of autologous platelet rich plasma (PRP) from a small sample of blood and for a preparation of a cell concentrate from bone marrow aspirate. They can be mixed with autograft and allograft bone prior to application to an orthopedic surgical site and can be used for wound care and other soft tissue applications as deemed necessary by the clinical use requirements.

WARNING AND PRECAUTIONS

1. Use proper safety precautions to guard against needle sticks.
2. Follow manufacturer instructions when using centrifuge. Use only EmCyte provided general purpose centrifuge. Outcomes using centrifuges from other manufacturers are unknown.
3. Do not use sterile components of this system if package is opened or damaged.
4. Single use device. Do not reuse.
5. Do not use after expiration date.
6. The surgeon is to be thoroughly familiar with the equipment and the surgical procedure prior to using this device.
7. Use prepared PPP or PRP within 4 hours after drawing blood from the patient, according to current AABB guidelines.
8. The patient is to be made aware of the general risks associated with drawing whole blood. These risks include, but are not limited to: hemorrhage, seroma formation, infection, and/or persistent pain at the site of aspiration.
9. The safety and effectiveness of this device for in vivo indications for use, such as bone healing and hemostasis, have not been established.
10. The PRP prepared by this device has not been evaluated for any clinical indications.
11. The PRP prepared by this device is NOT indicated for delivery to the patient's circulatory system.
12. Device is single use only. Do not attempt to clean or re-sterilize this product. After use, this product may be a potential biohazard.

POSSIBLE ADVERSE EFFECTS

1. Damage to blood vessels, hematoma, delayed wound healing and/or infection is associated with blood draw and/or surgical procedure.
2. Temporary or permanent nerve damage that may result in pain or numbness is associated with blood draw and/or surgical procedure.
3. Early or late postoperative infection is associated with surgical procedure.
4. Pain associated with site of whole blood harvest.

STERILITY

The platelet concentrating system kits are sterilized by ETO exposure. Do not use any component from an opened or damaged package. Do not resterilize. Single Use Only. Do not use past the expiration date.

INSTRUCTIONS FOR USE

NOTE: Use standard aseptic technique throughout the following procedure.

AbsolutePRP™ 20mL
GenesisCS Component Concentrating System

1. **WHOLE BLOOD DRAW:** Using the sterile filtered needle and 20mL syringe, draw 3mL of Citrate Anticoagulant into the 20mL syringe. Then slowly draw 17mL of blood from the patient filling the syringe to 20mL. Gently, but thoroughly mix the blood and anticoagulant upon collection to prevent coagulation.
2. **LOAD:** Unscrew the cap and slowly add the 20mL of anticoagulated whole blood into the Concentrating Device. Screw the cap onto the device and load into the Executive Series Centrifuge II machine or Platinum Series Machine.
3. **BALANCE:** Counterbalance the Concentrating Device. Place the counterbalance and Concentrating Device directly opposite to each other in the centrifuge rotor and make sure they contain the same amount of volume.
4. **SPIN:** Close centrifuge lid.
 - a. **For Executive Series Centrifuge:** Set RPM to 4000 and the time to 1 minute. Press the start button.
 - b. **For Platinum Series Centrifuge:** Toggle to the "AbsolutePRP Gold" programmed setting. Press the start button.
5. Once the centrifuge processing is complete, open centrifuge and remove Concentrating Device.
6. **PRP EXTRACTION:** After centrifugation, the platelets will remain suspended in the plasma solution. Using the 12mL syringe, aspirate the plasma portion containing the platelets, leaving the RBCs behind.

STORAGE CONDITIONS: Store in a cool and dry place.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Single
Use
Only

Authorized
Representative
Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands

EmCyte Corporation
4331 Veronica S.
Shoemaker Blvd., Suite 4
Fort Myers, Florida 33916
USA
Phone: 239-481-7725

Australian Sponsor
Emergo Australia
Level 20
Tower II, Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

EmCyte Corporation
4331 Veronica S. Shoemaker Blvd., Suite 4
Fort Myers, Florida 33916
USA

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


GS301-EU ** PLEASE DISCARD RED VENTED CAP FROM CONCENTRATING DEVICE BEFORE USE **

Note: Always swab self-sealing port with sterile alcohol prior to accessing with a sterile syringe.

<p>Step 1:</p>  <p>Using the sterile filtered needle, Draw 2mL of Citrate Anticoagulant into 20mL Syringe</p>	<p>Step 2:</p>  <p>Draw 18mL whole blood from the patient, filling the syringe to 20mL</p>	<p>Step 3:</p>  <p>Unscrew the cap and slowly add the 20mL of anticoagulated whole blood into the Concentrating Device</p>
<p>Step 4:</p> <p>Counterbalance the Concentrating Device according to the IFU insert and centrifuge at:</p> <p>For Executive Series Centrifuge</p> <p>1 minute 4000 RPM</p> <hr/> <p>For Platinum Series Centrifuge Toggle to AbsolutePRP™ Gold</p>	<p>Step 5:</p>  <p>High yielding PRP</p> <p>After centrifugation, the platelets will remain suspended in the plasma solution, forming the high yielding PRP.</p>	<p>Step 6:</p>  <p>Using the 12mL syringe, aspirate the high yielding PRP, leaving the RBCs behind.</p>