



HENSLER X[®] SURGICAL P R O D U C T S[™]

The choice is clear for this valuable bone source. The HBP was designed and tested to work as hard as the surgeons and support staff do; benefiting the patient and providing cost containment with every use.

HENSLER
BONE PRESS[™]



PRESS HEAD
COMPARISON




BONEBAC[®] PRESS

	Hensler Bone Press [™]	Bonebac [®] Press
Filtration/Collection	✓ Reinforced, proprietary, medical-grade stainless steel mesh	✗ No screen, porous holes into plastic disc
Filter surface area	✓ ~962 mm ²	✗ ~226 mm ²
Collection container size	✓ 80cc (x2)	✗ 30cc
Convenience	✓ Collection process staged out of surgical work area	✗ Hinders the operation between the scrub tech and surgeon
Ease of use	✓ Requires 5 actions	✗ Requires 10 actions
Press top	✓ Large, ergonomic top. Spring assisted return. Repeat press action easy & fast.	✗ Small top. Can be difficult to use. No spring return. Hard to repeat.


Kit contents:

- 1 Press top with spring return
- 1 Collection top
- 2 80 cc containers with ergonomic ribs
- 1 Angled tubing connector
- 1 3' 3/16" Sterile tubing
- Piece of mind



Kit contents:

- As shown below



DEMAND THE BEST.
HENSLER BONE PRESS[™]